

# 0 图文对齐训练

## 0.0 数据对齐

### 数据格式

论文中在这一阶段选择的数据集是 M3D-Cap。

链接：[3D多模态医疗数据集-图文对·数据集](#)

```
1  M3D_Cap/
2      ct_case/
3          000006/
4              Axial_non_contrast/
5                  0.jpeg
6                  1.jpeg
7                  .....
8                  text.txt
9      .....
10     ct_quizze/
11         000007/
12             Axial_non_contrast/
13                 0.png
14                 1.png
15                 .....
16                 text.txt
17     .....
18     .....
```

有点抽象的是里面的 3D 切片使用多个 jpeg 或 png 的 2D 切片组合起来的。

每个病例包含一个 **多个** 3D 切片和对应的文本描述，case 和 quizze 的区别是后者用于医学考试，质量更高。

文本示例：

```
1  Liver shows the following: cirrhotic changes with surface nodularity and hypertrophied caudate lobe multiple bi-lobar patches of low attenuation (<10 HU attenuation in the non contrast phase and <25 HU in the portal venous phase), suggestive of hepatic steatosis right hepatic lobe segment VI and left hepatic lobe segment IVa patches of contrast enhancement showing contrast enhancement at the early enhancement at the arterial and portal phases with no contrast wash out at the delayed phase, likely inflammatory in nature signs of fibrosis evident by surface nodularity, widened preportal and gallbladder fossae Portal hypertension and splenomegaly.Mild abdominal and pelvic ascites.
```

## 数据流

### 配置数据列表文件

示例:

```
1  {
2      "train": [
3          {
4              "image": "./data/M3D_Cap/ct_case/000001/Axial_non_contrast.nii.gz",
5              "text": "./data/M3D_Cap/ct_case/000001/report.txt"
6          },
7          {
8              "image": "./data/M3D_Cap/ct_case/000002/Axial_non_contrast.nii.gz",
9              "text": "./data/M3D_Cap/ct_case/000002/report.txt"
10         }
11     ],
12     "test": [
13         {
14             "image": "./data/M3D_Cap/ct_case/000003/Axial_non_contrast.nii.gz",
15             "text": "./data/M3D_Cap/ct_case/000003/report.txt"
16         }
17     ],
18     "val": [
19         {
20             "image": "./data/M3D_Cap/ct_case/000004/Axial_non_contrast.nii.gz",
21             "text": "./data/M3D_Cap/ct_case/000004/report.txt"
22         }
23     ]
24 }
```

### 数据预处理

图像数据增强及转换:

```
1  train_transform = mtf.Compose(                                # 数据增强
2      [
3          mtf.RandRotate90(prob=0.5, spatial_axes=(1, 2)),
4          mtf.RandFlip(prob=0.10, spatial_axis=0),
5          mtf.RandFlip(prob=0.10, spatial_axis=1),
6          mtf.RandFlip(prob=0.10, spatial_axis=2),
7          mtf.RandScaleIntensity(factors=0.1, prob=0.5),
8          mtf.RandShiftIntensity(offsets=0.1, prob=0.5),
9      ]
10     mtf.ToTensor(dtype=torch.float),                          # 数据类型转换
11 )
12
13
14 val_transform = mtf.Compose(                                   # 验证/测试 数据预处理
15     [
16         mtf.ToTensor(dtype=torch.float), # 数据类型转换
17     ]
18 )
19 set_track_meta(False)
```

文本数据截断：设置的 ClinicalBERT 的 max token 为 512，对于超出长度的文本，先以 '!' 为标志划分句子，然后在不超过长度上限的前提下随机挑选句子，组成新的长度限制内的文本描述。

原代码中的数据管道没有 resize 的操作，也就是默认所有输入的图片都是统一尺寸 (W, H, D) = (256, 256, 128)。因此，我在 CLIPDataset 类中新增了图片缩放的预处理，缩放结合了扫描的空间分辨率，确保处理后的扫描与物理空间相符。

## 返回格式

如下:

```
1  ret = {
2      'image': image,
3      'text': text,
4      'input_id': input_id,
5      'attention_mask': attention_mask,
6      'question_type': "Image_text_retrieval",
7  }
```

## OpenMind 数据适配

OpenMind 包含脑部 MRI 图像 3D 扫描，以及基本的脑部掩码等，文本数据只有简单的元数据（如性别、年龄、种族等）。为了进行图文训练，使用 Qwen3-VL-8B-Instruct 生成伪文本描述。

利用 Qwen3-VL 的视觉能力，我们输入 3D 切片的 3 个正交中心切面，以及相应的指令。示例如下：

**Visual Input (Montage: Axial | Coronal | Sagittal)**



输入 prompt:

```
1  Role: System
2  You are an expert Neuroradiologist. Your task is to provide a conservative, high-level
   radiological impression (summary) for the provided brain MRI case. Rules:
3  1. VOLUMETRIC SYNTHESIS: Although the input is a montage of representative orthogonal slices,
   interpret them as a complete 3D MRI volume.
4  2. INTEGRATE METADATA: Start the report by explicitly stating the scan modality and patient
   demographics.
5  3. PRIORITIZE MACRO-STRUCTURES: Assess symmetry, midline, and ventricular system.
6  4. BE CONSERVATIVE: Do not report minute findings unless they are undeniably abnormal.
7  5. TERMINOLOGY: Use 'Unremarkable' or 'Within normal limits'.
8  6. OUTPUT FORMAT: A single, flowing paragraph. No bullet points.
```

```

9
10 Role: User
11 [IMAGE_DATA_PLACEHOLDER] Analyze these representative orthogonal slices from a 3D volumetric
12 MRI scan. The views are presented as: Axial (left), Coronal (middle), Sagittal (right).
13 Case Metadata: Subject is a 26-year-old female. Modality: T1-weighted MRI. [Note: Image
14 quality is fair (Score: 2.3/5).] Note: Images are skull-stripped (background removed).
15 Clinical Indication: Routine follow-up for headache.
16 Please generate the Impressions.
17
18 Role: Assistant
19 Findings:
20

```

生成文本:


```

1 This 3D volumetric T1-weighted MRI of the brain in a 26-year-old female, despite suboptimal
2 image quality (score 2.3/5), demonstrates symmetric cerebral hemispheres with preserved sulcal
3 and gyral patterns across axial, coronal, and sagittal planes. The ventricular system appears
4 within normal limits in size and configuration, with no evidence of midline shift or
5 significant structural abnormality. No focal lesions, mass effect, or abnormal signal
6 intensity are identified. The overall impression is unremarkable.

```

## RadImageNet-alignment 数据适配

数据示例:



```

1 Image size: (224, 224), mode: RGB
2 Metadata: {'content_type': 'description', 'correct_text': None, 'is_abnormal': True,
3 'location': 'ankle foot', 'modality': 'mri', 'pathology': 'Plantar_plate_tear', 'question_id':
4 'description'}
5 [{'from': 'human', 'value': '<image>\nDescribe the medical imaging features and anatomical
6 details.'}]
7 {'from': 'template', 'value': 'magnetic resonance imaging of the ankle foot foot with plantar
8 plate tear present'}]

```

包含 metadata, 以及简单的描述。

原数据集共 750k 病例, 其中 MRI 600k 例, 覆盖不同部位与不同疾病, 数据总览如下:

```

1 === Modality Counts Before Cleaning ===
2 Modalities Total: {'ct': 228902, 'mri': 604775}
3 === Modality Counts After Cleaning ===

```

```

4  Modalities Total After: {'mri': 604775}
5
6  === Overall Metadata Statistics ===
7  Field: content_type
8      description: 604775
9
10 Field: correct_text
11     None: 604775
12
13 Field: is_abnormal
14     True: 482222
15     False: 122553
16
17 Field: location
18     ankle foot: 162517
19     brain: 40001
20     hip: 46477
21     knee: 161744
22     abdomen: 82563
23     shoulder: 47206
24     spine: 64267
25
26 Field: modality
27     mri: 604775
28
29 Field: pathology
30     Plantar_plate_tear: 599
31     achilles_pathology_: 8277
32     atfl_pathology: 10819
33     bone_inflammation: 36099
34     cfl_pathology: 3851
35     chondral_abnormality: 85073
36     coalition: 58
37     deltoid_pathology: 2478
38     extensor_pathology_: 618
39     fat_containing_tumor: 244
40     flexor_pathology_: 3259
41     hematoma: 747
42     intra: 5954
43     lisfranc_pathology: 47
44     normal: 122553
45     osseous_disruption: 11826
46     osseous_neoplasm: 1317
47     peroneal_pathology: 6589
48     plantar_fascia_pathology: 5752
49     post_op: 5659
50     soft_tissue_edema: 15691
51     soft_tissue_fluid: 37215
52     soft_tissue_mass: 8472
53     spring_ligament_injury: 97
54     syndesmosis_pathology: 343
55     acute_infarct: 448
56     arteriovenous_anomaly: 241

```

57 chronic\_infarct: 2040  
58 edema: 111  
59 extra: 1113  
60 focal\_flair\_hyper: 672  
61 pituitary\_lesion: 69  
62 white\_matter\_changes: 9185  
63 abductor\_pathology\_: 765  
64 capsular\_pathology: 109  
65 chondral\_pathology: 7260  
66 hamstring\_pathology: 212  
67 labral\_pathology: 36163  
68 marrow\_inflammation: 4295  
69 osseous\_lesion: 2364  
70 acl\_pathology: 9059  
71 fcl\_pathology: 408  
72 fracture: 1402  
73 mcl\_pathology: 7123  
74 meniscal\_abnormality: 40498  
75 muscle\_strain: 464  
76 patella\_pathology: 1266  
77 pcl\_pathology: 667  
78 post\_operative\_acl: 77  
79 quad\_pathology: 450  
80 soft\_tissue\_fluid\_collection: 23060  
81 adrenal\_pathology: 581  
82 arterial\_pathology: 76  
83 ascites: 699  
84 bil\_dil: 216  
85 bladder\_pathology: 130  
86 bowel\_abnormality: 151  
87 bowel\_inflammation: 239  
88 bowel\_mass: 93  
89 degenerative\_changes: 165  
90 dilated\_urinary\_tract: 55  
91 enlarged\_organ: 183  
92 gallbladder\_pathology: 434  
93 intraperitoneal\_mass: 804  
94 liver\_disease\_: 971  
95 liver\_lesion: 2987  
96 marrow\_abn: 649  
97 ovarian\_pathology: 1673  
98 pancreatic\_lesion: 996  
99 prostate\_lesion: 65  
100 renal\_lesion: 1834  
101 soft\_tissue\_collection: 161  
102 splenic\_lesion: 214  
103 uterine\_pathology: 3219  
104 acj\_oa: 2392  
105 biceps\_pathology: 1450  
106 ca++\_tendinosis: 69  
107 ghj\_oa: 2618  
108 infraspinatus\_pathology: 80  
109 subscapularis\_pathology: 128

```

110     supraspinatus_pathology: 4680
111     cord_pathology_: 642
112     cystic_lesions: 1221
113     disc_pathology: 29370
114     dural_epidural_abn: 7410
115     facet_arthropathy: 325
116     foraminal_pathlogy: 10735
117     osseous_abn: 2095
118     scoliosis: 1807
119
120     Field: question_id
121     description: 604775

```

对于 2D 图片，将其复制 32 层置于中心，其余全黑，**在平均池化的时候乘上缩放因子**  $128/32 = 4$ ，避免信号过弱。32 层是因为模型最大核大小为 7，深度下采样 16 倍，选择 16 或 32 可以较好保持信息。

对于文本，采用 metadata 和原文本结合的方式构造。

## 训练

### 训练配置

```

1  #!/bin/bash
2  deepspeed src/train/train_clip.py \
3      --deepspeed ./scripts/zero2.json \
4      --language_model_name_or_path medicalai/ClinicalBERT \
5      --wb_name DCFormer_SigLIP \
6      --vision_encoder "dcformer" \
7      --loss_type "sigmoid" \
8      --data_root ./data \
9      --cap_data_path ./data/datalist.json \
10     --max_length 512 \
11     --bf16 True \
12     --output_dir ./output/DCFormer_SigLIP \
13     --num_train_epochs 100 \
14     --per_device_train_batch_size 42 \
15     --per_device_eval_batch_size 4 \
16     --gradient_accumulation_steps 1 \
17     --eval_strategy "no" \
18     --eval_accumulation_steps 1 \
19     --eval_steps 0.04 \
20     --save_strategy "steps" \
21     --save_steps 1000 \
22     --save_total_limit 1 \
23     --learning_rate 1e-4 \
24     --weight_decay 0.1 \
25     --warmup_ratio 0.03 \
26     --lr_scheduler_type "cosine" \
27     --logging_steps 0.001 \
28     --gradient_checkpointing False \
29     --dataloader_pin_memory False \

```

唯一更改的地方: `per_device_train_batch_size`: 64 -> 42

## 训练进度

| 数据集                   | 训练(cases)   | 图像类型 | 文本  |
|-----------------------|-------------|------|-----|
| OpenMind              | 10000 (10%) | 3D   | 伪文本 |
| RadImageNet-alignment | 12000 (2%)  | 2D   | 原生  |
|                       |             |      |     |

训练监控: [MRI-3D-VLM Table – Weights & Biases](#)